Correlation of Thin and Bulk Hydrogel Mechanics to Measure the Mechanical Properties of Soft

Biomaterials

By: Nathan Richman

Advisor: Professor Melih Eriten

Department of Mechanical Engineering, University of Wisconsin-Madison

Abstract


Bacterial biofilms are complex structures of bacteria that can form on nearly every exposed surface¹. While biofilms negatively impact many aspects of industry and are associated with nearly 80% of infections in humans, humans wouldn't be able to survive without them^{2,3}. Understanding the prevention and growth of biofilms has been a topic in research for years. Recent research has begun to quantify mechanical properties of biofilms and their responses to mechanical stress, however there are limited techniques available to accurately characterize these properties⁴. This study investigates macro and micro scale characterization of thin and bulk gel layers using gelatin, PDMS and polyacrylamide to model the mechanical properties of bacterial biofilms. The methods and results of this study can be used to characterize the mechanical properties of bacterial biofilms grown in a variety of conditions.

Background and Motivation

Hydrogel Mechanics

Poroelasticity:

- Time and area dependent stress for a given displacement due to diffusion (**D**) limited migration of incompressible solvent⁵⁻⁷
- As a result, much higher initial force F₀ vs. equilibrium Force F_m



Viscoelasticity:

- Time dependent reduction in stress for a constant displacement due to molecular movement⁷
- Polymer chains "store" elastic energy, but the movement of molecules leads to energy being "lost" to heat.
- E'/G' → Storage Modulus

of substrate effects:

- E"/G" → Loss Modulus
- $G^* = G' + iG''$

Problem:

Hard to test mechanical properties of thin films accurately because

Mechanical testing of bulk specimen

G*, v

Figure 2: Illustration of chemistry and molecular physics behind viscoelasticity

(Left): Resting state of polymer chains intertwined. (Center): Immediately after

force is applied the chains deform somewhat. (Right): After time the chains

 $\oint F_0 \approx \frac{16}{2} R^{1/2} d^{3/2} G \quad (1)^5$

____ d[

Diffusion limited

migration of solvent

Indenter

←---

become detangled and relax some stress

2-15mm

Much larger sample

compression

Larger forces (1-100N)

Minimal substrate effects

Indentation, rheometer, bulk

Figure 1: Basic formula

nanoindentation with

Indenter of radius R

indentation depth d

and gel height h

and properties of

Mechanical testing of thin film G*, v, D_w h=30-200μm E^1 , v^1

- Typical size of a biofilm
- Very low forces (μN-mN)
- Substrate effects
- Nanoindentation and AFM

$$F_0 \approx 16/3 * R^{1/2} \delta^{3/2} G * f(\sqrt{Rh}/d)$$
 (2)⁸

$$F_{\infty} \approx \frac{8R^{1/2}\delta^{3/2}G}{3(1-v)} * f(\sqrt{Rh}/d)$$
 (3)⁸

Hypothesis:

We can develop methods for testing mechanical properties of a thin hydrogel and verify results based on properties of bulk hydrogel. These results will provide a verified method for determining mechanical properties of soft biomaterials

Results

Bulk Gel (Gelatin):




Figure 5: Stress plotted against stretch showing the strain hardening properties of gelatin at various concentrations. All data is fit with the 1st order Ogden model. Equation for fit is shown on the graph.

Table 1: Average μ₁ and α₁ for Gelatin of 4.68 (.52) 2.46 (.04) 11.51 (.52) varying concentrations. General trend seemed 4.21 (3.64) 2.83 (.45) 11.91 (3.6) to be decreasing μ .67 (.42) 3.37 (.58) 2.25 (.71) Standard deviation in

Poro-viscoelasticity: Minute scale





Figure 6: Normalized Force plotted against Normalized time (0-5min) Separation of the curves at this time scale indicates that the viscoelasti properties of the matrix dominated over the poroelastic properties.

Nanoindentation (Gelatin):






Figure 7: Fitted force relaxation for a poroelastic gel, however the gel was not thin. The data was fit to two models, either the model presented in introduction⁸, or a bulk poroelastic model⁶. The function was plotted, and Dt/Rh was calculated using various values of D until the fit was sufficient. Possible fractures indicated

Displacement Rate Dependence:




Figure 8: Force vs. Depth at two strain rates. Again samples from different locations (same color) showed very similar mechanical properties. Secondly, the increased apparent stiffness at the faster strain rate (blue) indicates poro-viscoelastic properties, but it is not possible to decouple them with this method.

Discussion of Results

Presentation of Property	Bulk Gel	Nanoindentation	Explanation
Hyperelasticity			 Difficult to identify hyperelasticity or the nanoscale due to fractures, and indenter and instrument limits.
Poro- viscoelasticity			 Expect to see poro-viscoelasticity at any length scale Time scale length vs. poroelasticity: is poroelasticity present? Strain rate dependence observed as expected.

Poroelasticity:

- Poroelastic model fits fairly well, but still unsure whether poroelasticity is dominant.
- Diffusion very small, possibly supports poroelasticity
- Need to do stress relaxation at greater time scale (hours)

Bulk vs. Nanoscale:

- Shear modulus on same order of magnitude on nano and macro scale (5.83kPa, n=1, vs. 9.98kPa, n=2
- Viscoelasticity grossly seen on both length scales
- Hyperelasticity not seen, maybe too small of length scale to be seen




Figure 9: Illustration showing how poro and viscoelasticity can be separated. Poroelastic effects show normalized force vs. time as distinct curves, but normalizing time to the length scale collapses the curves into one. Adapted from Wang et al. 2014¹⁰

Limitations:

- Possible fracture due to drying: Gelatin hydrogel has water evaporate from the surface and so the surface tends to be more brittle than the bulk.
- Not thin indentation: thick (1-2mm) gel was used for indentation.
- Greater implications of drying: Surface characteristics completely different than bulk possibly


Conclusion and Future Work

Conclusions:

- Shear moduli fairly close at each scale, non-quantified properties qualitatively similar. More work needed to draw strong conclusions.
- Hard to conclude about similarity of properties due to drying and not using thin indentation.

Future Work:

- Replicate Diffusivity at different values of \sqrt{Rh}/d , at different values, should get same diffusivity:
- Use crosslinked gel, immerse in solvent
- Eventually use to measure bulk properties of biomaterials



Acknowledgements

- Thank you to the Materials Science Center (MSC) for the use of the Hysitron Triboindenter
- Thank you to the Biomaterials Teaching Lab for the use of the MTS machine and rheometer
- Thank you to Professor Melih Eriten, Guebum Han, and Dr. Utku Boz for their continued guidance throughout the duration of the project and for the initial project idea.

References

- 1. D. Pavithra and M. Doble, *Biomed. Mater.*, vol. 3, no. 3, p. 34003-16, 2008. 2. K. Epstein, I. Hochbaum, P. Kim, and J. Aizenberg, Nanotechnology, vol. 22, no. 49, p. 494007 (8pp), 2011.
- 3. C. H. Thomson, Int. Wound J., vol. 8, no. 1, pp. 63-67, 2011. 4. D. Head, Adv. Exp. Med. Biol, vol. 915, pp. 99-114, 2016.
- 5. E. P. Chan, Y. Hu, P. Johnson, Z. Suo, C. M. Stafford, *Soft Mat.*, vol. 8, no. 5, p. 1492, 2012. 6. Y. Hu, X. Zhao, J. J. Vlassak, and Z. Suo, *Appl. Phys. Lett.*, vol. 96, no. 12, 2010.
- 7. T. Osswald, G. Menges, Materials Science of Polymers for Engineers, Hanser Publishers, 2003.
- 8. Y. Hu, E. P. Chan, J. J. Vlassak, and Z. Suo, J. Appl. Phys., vol. 110, no. 8, 2011. 9. R. W. Ogden, Proceedings of the Royal Society of London, vol. 326, no. 1567, pp. 565–584, 1972. 10. Q. M. Wang, A. C. Mohan, M. L. Oyen, and X. H. Zhao, *Acta Mech. Sin. Xuebao*, vol. 30, no. 1, pp. 20–27, 2014.

Bulk Gel:

Methods

- Varying stiffness gels cast into 14-50mm diameter molds
- Compressed and rotationally sheared between two parallel plates.





Figure 4 (above): Setup of bulk gel testing, the plates can either move along or rotate about the z-axis, 10mm Scale bar. Figure 5 (right): Image of surface of gel after nanoindentations, these spots show possible fractures. Scale bar shown is 50µm.

 $G' = \frac{\tau_0}{\nu_0} \cos(\delta) (5)^7$

Thin Film:

- . Coat petri dish with thin gel
- 2. Indent 10-30μm, record force and time 3. Define τ and x to be dimensionless time and
- displacement: $\tau \equiv \frac{1}{Rh}$, $x \equiv \frac{1}{h}$
- 4. Plot dimensionless force with respect to these variables and use curve fitting with n = 2 to find α , β, and γ.







Figure 3: Properties and general testing strategy in thin film and bulk material. (Left) From indentation and shear, one can get the shear modulus, Poisson's ratio and diffusivity of water. (Right) From compression and rotational shear one can get shear modulus, and Poisson's Ratio.

parentheses, n=2.

- Stress/strain: can find G through eq. 4



 $G'' = \frac{\tau_0}{100} \cos(\delta) (6)^7$